

# CBSPD Practice Test Questions and Answers

## 1. What is the purpose of leak testing a flexible endoscope?

- A) To verify sterilization efficacy
- B) To detect internal damage that could allow fluid ingress
- C) To measure channel diameter
- D) To confirm cleaning chemical compatibility

## 2. What is an Automated Endoscope Reprocessor (AER) used for?

- A) Sterilizing rigid scopes via steam
- B) Automating the HLD rinse and drying cycles for flexible endoscopes
- C) Performing only leak testing
- D) Replacing manual cleaning steps entirely

## 3. How should flexible endoscopes be stored after reprocessing?

- A) Coiled tightly in a sealed plastic bag
- B) Hanging vertically in a ventilated cabinet with caps removed
- C) Submerged in disinfectant solution
- D) Wrapped in sterile drapes in a drawer

## 4. When should the integrity of the sterile packaging be inspected?

- A) Only when there is visible damage
- B) Before and after sterilization, and before use
- C) Only after sterilization
- D) Once a month

**Answers: 1-B 2-B 3-B 4-B**

For More CBSPD Questions and Answers FREE, CBSPD Online Prep Training, CBSPD Exam, CBSPD Study Guide, CBSPD Flashcards, CBSPD Quizzes visit:

## CBSPD Practice Test