

CSSD Practice Test Questions and Answers

1. What does CSSD stand for in healthcare sterilization?

- A) Central Sterile Supply Department
- B) Certified Surgical Service Dispatcher
- C) Clinical Sterilization Standards Division
- D) Certified Sterile Service Distributor

2. What is the recommended temperature and time for steam sterilization of wrapped instruments?

- A) 250°F for 15 minutes
- B) 270°F for 30 minutes
- C) 275°F for 15 minutes
- D) 285°F for 10 minutes

3. Which method is most effective for sterilizing heat-sensitive medical devices?

- A) Steam sterilization
- B) Ethylene oxide (EtO) sterilization
- C) Dry heat sterilization
- D) Boiling water

4. What is the proper sequence for processing reusable medical instruments?

- A) Sterilization, cleaning, packaging, storage
- B) Packaging, cleaning, sterilization, storage
- C) Cleaning, inspection, packaging, sterilization, storage
- D) Storage, cleaning, sterilization, packaging

Answers: 1-A 2-C 3-B 4-C

For More CSSD Questions and Answers FREE, CSSD Online Prep Training, CSSD Exam, CSSD Study Guide, CSSD Flashcards, CSSD Quizzes visit:

CSSD Practice Test