

# BE/B Eng Bachelor of Engineering Practice Test Questions and Answers

## 1. EKG is short for

- A) Electrtokinetcography
- B) Electromyography
- C) Electroencephalography
- D) Electrocardiography

## 2. Paste for electrodes

- A) absorbs contact impedance
- B) reduces contact impedance
- C) equates contact impedance
- D) increases contact impedance

## 3. Muscle palsies indicate

- A) involuntary muscles
- B) voluntary muscles
- C) active muscles
- D) paralyzed muscles

## 4. Buffer amplifier transforms

- A) high impedance signals to low impedance signals
- B) low impedance signals to high impedance signals
- C) dc impedance signals to ac impedance signals
- D) ac impedance signals to dc impedance signals

**Answers: 1-D 2-B 3-D 4-A**

For More BE/B Eng Bachelor of Engineering Questions and Answers FREE, BE/B Eng Bachelor of Engineering Online Prep Training, BE/B Eng Bachelor of Engineering Exam, BE/B Eng Bachelor of Engineering Study Guide, BE/B Eng Bachelor of Engineering Flashcards, BE/B Eng Bachelor of Engineering Quizzes visit:

## BE/B Eng Bachelor of Engineering Practice Test